

SOFT CONTACT LENSES: INSTRUCTION FOR USE

A Guide to Soft Contact Lens Wear



Fig. 1



Fig. 2



Fig. 3



Fig. 4

Sponsor:

CooperVision Australia Pty Ltd
115 Hayward Avenue, Torrensville
SA 5031, Australia

Legal Manufacturer:


CooperVision Manufacturing Limited,
South Point, Hamble, Southampton,
SO31 4RF, UK



Australia: www.coopervision.net.au
New Zealand: www.coopervision.co.nz

Soft contact lenses in buffered saline.

<p>Indications and Intended Use / Purpose: For the correction of ametropia. Clinical Benefit and Performance Characteristic: Improving visual acuity while wearing the contact lenses. Toric lenses correct astigmatism. Multifocal lenses correct presbyopia. Target User Group / Intended Users: For users with ametropia (hyperopia, myopia, astigmatism, presbyopia) with healthy eyes. There is no restriction on age; suitability should be determined by an eye care professional.</p>	
Material Family Name Description	Wear Schedule
<p>Frequent Replacement Contact Lenses (Not Daily Disposable) Manufacturer's Recommended Replacement Schedule: One Month</p>	
<p>comfilcon A: Sphere/XR; Asphere; Multifocal; Toric; XR Toric; Toric Multifocal</p>	<p>Daily Wear or Extended Wear Suitability should be determined by an eye care professional.</p>
<p>fanfilcon A: Sphere; Toric</p>	
<p>somofilcon A: Asphere; Multifocal; Toric</p>	
<p>omafilcon B: Sphere; Multifocal; Multifocal XR; Multifocal Toric; Toric; Toric XR</p>	<p>Daily Wear</p>
<p>Daily Disposable Contact Lenses Manufacturer's Recommended Replacement Schedule: Daily Disposable</p>	
<p>ocufilcon D: Sphere</p>	<p>Daily Wear</p>
<p>somofilcon A: Sphere; Asphere; Multifocal; Toric</p>	
<p>stenfilcon A: Asphere; Multifocal; Toric</p>	
<p>omafilcon A: Sphere; Multifocal</p>	
<p>riofilcon A: Sphere</p>	
<p>Indications: For the control of myopia. Intended Use / Purpose: For the correction and control of myopia. Intended Clinical Benefits and Performance Characteristic: Improving visual acuity and slowing the progression of myopia while wearing contact lenses. Target User Group / Intended Users: For children with healthy eyes within the age range of 6-18 years, who require myopia correction and want to slow the progression of myopia. Suitability should be determined by an eye care professional.</p>	
Material Family Name Description	Wear Schedule
<p>Daily Disposable Contact Lens</p>	
<p>omafilcon A: Dual Focus</p>	<p>Daily Wear</p>

 Do not use if the product sterile barrier system or its packaging is compromised. The sterile barrier for contact lenses is the intact blister seal. Always wash and dry your hands before handling your contact lenses.

Lens Application:

1. Before application, ensure the lens is not inside out. Place on your index finger to verify the shape (Fig. 1 & 2).
2. Using the middle finger, hold your lower lid down. With the index finger of the other hand, raise upper lid. Place the lens onto the eye (Fig. 3). Repeat for the other eye.

Lens Removal:

Separate the lids with the middle finger of each hand. Look up and with the index finger of the lower hand; slide the lens down onto the white of the eye. Pinch the lens and lift out with the thumb and index finger (Fig. 4). Repeat for the other eye.

Disinfection/Soaking for Frequent Replacement Contact Lenses:


For lenses worn more than once (single wearer, multiple use) clean and disinfect them strictly following the instructions on your lens care solution bottle or as directed by your eye care professional (for example: ophthalmologist, optometrist, or optician). Follow instructions on the lens care products for proper contact lens case care, including rinsing the lens case with fresh solution and wiping it with a clean tissue or air drying after use. Cases should be replaced as advised on the lens care product instructions.

If Contact Lenses Dry Out:

If contact lenses need rewetting during wear, certain lens care products (used with frequent replacement lenses as indicated) can be used to rewet lenses, or lubricating/rewetting/comfort drops used with soft contact lenses, as recommended by your eye care professional. However, if a lens becomes completely dried out, it should be discarded and replaced with a new one.

Daily Disposable / Single Use Contact Lenses:

Daily disposable lenses are intended to be worn only once, for daily wear, and then discarded. They should not be worn while sleeping. Daily disposable lenses are not designed to be cleaned or disinfected. Daily disposable lenses worn more than once may become damaged or contaminated and could put you at risk of eye infections, inflammation, discomfort, or other problems.

 Symbol on labelling indicates the contact lens is for daily disposable wear and not to reuse. When symbol is not present, the lens is indicated for other wear schedule.

Precautions:

1. For all types of contact lenses, your eye care professional will recommend lens selection, wearing, cleaning and replacement regime. Do not vary these. If in doubt seek the advice of your eye care professional.
2. Contact lenses are medical devices and should always be fitted by your eye care professional.
3. If your eyes become red or irritated, or you experience any pain or unexpected change in vision, remove your lenses immediately and consult your eye care professional.
4. Avoid lens wear if you are unwell, in particular with any cold or flu-like symptoms.
5. Avoid contaminating lenses with make-up, aerosol sprays or wearing them near noxious fumes.
6. Contact lenses should be removed immediately after contact with noxious vapour, e.g. chemical or hazardous substances, or hazardous environment with ocular impact.
7. Replace damaged, contaminated or scratched lenses immediately.
8. Avoid long nails, which harbour microbes and could damage the eye or lens. A scratched eye is more vulnerable to infection.
9. Keep lenses out of the reach of children who are not users.
10. Exposure with water (e.g. from the tap, bottled, bath, pool, ocean or during water sports) increases the risk of microbial infection, so do not expose your contact lenses to water, including while wearing them. If they are exposed, the lenses should be discarded and replaced with a new pair.
11. Do not change from the lens care solutions recommended by your eye care professional.
12. Do not change lens type or parameters without consulting your eye care professional.
13. Always have a pair of glasses with a current prescription available to use if you become unable to wear your contact lenses, or where contact lens wear is not advised.

Additional precaution for Frequent Replacement Lenses

14. For lenses worn more than once, always clean and disinfect your lenses as instructed, after lens removal.
15. If wearing lenses on an extended wear basis, remove and dispose of your lenses the evening before you are due to replace them, inserting new lenses the following day.

Warnings:

1. **DO NOT** sleep in your lenses unless your eye care professional has advised it is safe to do so.
2. **DO NOT** wear your lenses beyond the period recommended by your eye care professional.
3. **DO NOT** use household products (e.g. disinfectants) on your lenses.
4. **DO NOT** wear your lenses when using eye drops (unless advised to do so by a medical professional).
5. **DO NOT** use saliva to wet your lenses.
6. **DO NOT** expose your contact lenses to water while you are wearing them.
7. **DO NOT** rinse or handle your lenses or lens case with or expose them to water as it increases the risk of serious microbial infection.
8. **DO NOT** use if the expiry date has passed.
9. **DO NOT** share your contact lenses. Sharing lenses greatly increases the chance of eye infections.

Contraindications (reasons not to use)

Do not use contact lenses when any of the following conditions exist:

- Acute and subacute inflammation or infection of the anterior chamber of the eye.
- Any eye disease, injury, or abnormality that affects the cornea, conjunctiva, or eyelids.
- Corneal hypoesthesia (reduced corneal sensitivity), if not aphakic (natural lens of the eye is present).
- Severe insufficiency of lacrimal secretion (dry eyes).
- Any systemic disease that may affect the eye or be exaggerated by wearing contact lenses.
- Allergic reactions of ocular surfaces or adnexa (eyelids and lashes) that may be induced or exaggerated by wearing contact lenses or use of contact lens solutions.
- Any active corneal infection (bacterial, fungal, protozoal or viral).
- If eyes become red or irritated.
- Use of a medication that is contraindicated, including medications for the eye.

Adverse Reactions (Undesirable Side Effects)

You may become aware of the following symptoms while wearing contact lenses:

- Reduced lens comfort or feeling of something in the eye
- Burning, grittiness, itchiness, dryness of the eyes
- Watery eyes and/or unusual eye discharge
- Sensitivity to light (photophobia)
- Swelling or inflammation of the eye (swollen eyelid, eye redness)
- Temporary decrease in vision

What to do when symptoms occur: You should immediately remove the lenses. If the symptoms persist, a serious eye condition may be present. Seek prompt care from your eye care professional so that any potential medical condition can be identified and treated, if necessary, in order to avoid serious eye damage.














Medical conditions that may be associated with these symptoms include: corneal abrasion, corneal oedema (swelling), infiltrative keratitis (inflammation of the cornea), corneal infection or ulcer, corneal neovascularisation, and papillary conjunctivitis (eyelid inflammation).

Disposal: use standard waste, recycle whenever possible.

Any serious incident that has occurred in relation to the product should be reported to the manufacturer and applicable competent authority.

SYMBOLS KEY:

The following symbols may appear on the label or carton.

SYMBOL	DEFINITION
	Caution: By law, this device must be prescribed by a registered practitioner
	Caution / See Instructions for Wearers
	Use by Date (expiration date)
	Batch Code
	Sterile using Steam Heat
	UV Blocking (if applicable)
	Manufacturer
	Authorized representative in the European Community
	Medical device
	Do not use if package is damaged
	Consult instructions for use
	Do not re-use (if applicable)
	Date of Manufacture